

KNOW YOUR HEALTH.

EMPOWER YOUR LIFE.



BIO-VITALS KIOSK

Smart Health Screening
Solution for instant health insights.

Model: KH-0080



Height | Weight Measurement



BMI Analysis with 18 parameters



Blood Pressure & Body Temperature



Blood Oxygen Level (SpO₂)



Heart Rate Tracking



Smart Interactive Display



Instant Health Reports

IDEAL FOR



Clinics/
OPDS



Metro / Bus
Stations



Gyms



Corporate/Offices
Factories



Pharmacies



Wellness
Centers



Schools/
Colleges



Government
Institutions



Housing
Societies



Recruitment
Camps

FAST ACCURATE INSTANT

A SMART STEP TOWARDS BETTER HEALTH

Monitor essential health parameters instantly with **ONYX Health+**. Designed for quick, accurate, and user-friendly preventive health screening.


WHY ONYX Health+

- ✓ Promoting Preventive Health Care
- ✓ Supporting Community Health
- ✓ Quick Results
- ✓ Modern Design
- ✓ Data in Excel QR code format

TECHNICAL SPECIFICATIONS

Parameter	Specification
Model	KH-0080
Product Name	Bio-Vitals kiosk
Height Measurement Range	80 cm – 210 cm (Accuracy: ±0.5 cm)
Weight Capacity	Up to 180 kg (Accuracy: ±100 g)
BMI Calculation	Automatic
BP Measurement Method	Oscillometric (Accuracy: ±3 mmHg)
Display	10.1" LCD Touch Screen
Language	English
Connectivity	WiFi Module (Optional SIM)
Power Supply	220V AC / 50Hz
Power Consumption	10W
UPS (Optional)	600 VA
Printer	Thermal Printer
Body Material	Mild Steel (MS) with Powder Coated Finish
Operating & Storage Temperature	10°C to 40°C
Dimensions (W x D x H) (approx.)	Body: 300 x 85 x 2100 mm Base: 500 x 25 x 800 mm
Net Weight	49 Kgs.
Gross Weight	53 Kgs

Specifications are subject to change without notice.



*kiosk print-out shows all other BMI parameters also

How It Works

1. Stand on the Machine
2. Checks Height, Weight, Blood Pressure (BP), SpO₂ & Body Temperature
3. View Your Health Report on the Screen
4. Get an Instant Printout
5. Scan & Download Your Report via QR Code (Optional)
6. Download Bulk Reports in Excel Format Through the Portal
7. API Integration Available
8. UPI QR Code Payment Mode (Optional)

📍 ONYX EQUIPMENTS (INDIA) PVT.LTD. - PUNE, DELHI, KOCHI.

☎ 81491 22926/ 82370 06614

✉ sales@onyxhealthplus.com

🌐 www.onyxhealthplus.com