

KNOW YOUR HEALTH.

EMPOWER YOUR LIFE.



BMI - BP ANALYZER

Smart self-service health screening solution for instant health insights.

Model: KH-0060

 Height | Weight Measurement

 BMI Analysis with 18 parameters

 Blood Pressure Monitoring

 Heart Rate Tracking

 Smart Interactive Display

 Instant Health Reports

IDEAL FOR

 Clinics / Hospitals

 Metro / Bus Stations

 Gyms

 Corporate Offices

 Pharmacies

 Wellness Centers

 Schools / Colleges

 Government Institutions


 **FAST.** |  **ACCURATE.** |  **EASY.**

A SMART STEP TOWARDS BETTER HEALTH

Monitor essential health parameters instantly with **ONYX Health+**. Designed for quick, accurate, and user-friendly preventive health screening.

WHY ONYX Health+

- ✓ Promoting Preventive Health Care
- ✓ Supporting Community Health
- ✓ Quick Results
- ✓ Modern Design
- ✓ Affordable



TECHNICAL SPECIFICATIONS

Parameter	Specification
Model	KH-0060
Product Name	BMI-BP Analyzer
Height Measurement Range	80 cm – 210 cm (Accuracy: ± 0.5 cm)
Weight Capacity	Upto 180 kg (Accuracy: ± 100 g)
BMI Calculation	Automatic
BP Measurement Method	Oscillometric (Accuracy: ± 3 mmHg)
Display	10.1" LCD Touch Screen
Language	English
Body Material	Mild Steel (MS) with Powder Coated Finish
Power Supply	220V AC / 50Hz
Power Consumption	10W
UPS (Optional)	600 VA
Printer	Thermal Printer
Connectivity	SIM module (Optional Wi-Fi)
Data Storage	Not available
Operating & Storage Temperature	10°C to 40°C
Dimensions (W x D x H) (approx.)	Body: 300 x 85 x 2100 mm Base: 500 x 25 x 800 mm
Net Weight	49 Kgs.
Gross Weight	53 Kgs

Specifications are subject to change without notice.

HOW IT WORKS

